

Skin Sensitization Study in Guinea Pigs (Maximization Test)

Master Laboratory Co., Ltd. Animal laboratory

Address: 3F, No. 221, Sec. 1, Zhongxing Rd., Zhudong Township,

Hsinchu County 31053, Taiwan (R.O.C.) TEL: (886-2)25176117



Report No.: MSA-202010-409-T03

Ear-loop Skin Sensitization Study in Guinea Pigs (Maximization Test) STUDY REPORT

Sponsor: CHIN HSIUNG FIBER CO., LTD

Testing Institution: Master Laboratory Co., Ltd.

November 2020



Report No.: MSA-202010-409-T03

Report No.: MSA-202010-409-T03

Test article registration date: 10.22.2020

Experimental starting date: 10.27.2020

Animal in-housing: 10.27.2020

Extraction of test article date: 10.30.2020, 11.06.2020, 11.20.2020

Test article administration during induction phase: 11.02.2020, 11.09.2020

Test article administration during challenge phase: 11.23.2020

Observation of skin response: 11.25.2020-11.26.2020

Animal sacrifice date: 11.26.2020

Experimental complete date: 11.26.2020

Study Announcement

- 1. The study report is valid for the test article used only, and shall not be partly recopied or extracted for another object.
- 2. The study report is invalid without the endorsement of Master Laboratory Co., Ltd.

Report No.: MSA-202010-409-T03

SIGNATURE OF STUDY PERSONNEL

Study Director

130 Man Meang

Facility Management



Report No.: MSA-202010-409-T03

GLP COMPLIANCE STATEMENT

The study met with the technical requirements of the protocol, and all applicable guidance and regulations, which included the Good Laboratory Practice for Non-clinical Laboratory Studies (FDA, 21 CFR, Part 58, 2019) and Good Laboratory Practice for Non-clinical Laboratory Studies (Food and Drug Administration, R.O.C., 2019). There were no deviations from the approved study plan and no adverse problem that would affect the integrity of this study or the interpretation of the study result. Because the test article is a proprietary product of the sponsor, all the contents related with test article in 21 CFR Part 58 (US FDA) are not applicable to this study (21 CFR §58.105, §58.113, FDA).

Study Director

Bo Han Huang

Date



Inspection record:

FINAL REPORT

Report No.: MSA-202010-409-T03

QUALITY ASSURANCE STATEMENT

To comply with the "Good Laboratory Practice for Nonclinical Laboratory Study", Quality Assurance Department has audited the facility, equipment, personnel, test methods, raw data, and records regularly.

The study report has been reviewed and approved. The experiments were conducted according to the protocol. All original records, raw data, and documents are truthfully transferred and addressed in the results of this report.

Inspection Contents	Date of inspection
Before the test (test execution protocol, requisitions, contracts)	10.27.2020
Test (test substance data sheet, animal quarantine, standard operating procedures)	11.02.2020

After the test (complete the original data, report reviews).....

Quality Assurance unit in charge

Ying Chun Chen Ying Chun Chen

Date

11.27.2020



Report No.: MSA-202010-409-T03

Contract Research Organization and Sponsor Information

1. CRO:

a. Title:

Master Laboratory Co., Ltd.-Animal laboratory

b. Address:

3F, No. 221, Sec. 1, Zhongxing Rd., Zhudong Township, Hsinchu

County 31053, Taiwan, R.O.C.

2. Sponsor:

a. Title:

CHIN HSIUNG FIBER CO., LTD

b. Address:

NO.37-12, Ta Hsin RD., Pu Yen Hsiang, Chang Hua Hsien, Taiwan



Report No.: MSA-202010-409-T03

CONTENTS

Summary	8
Introduction	9
Materials and Methods	
Results	17
Conclusion	17
References	18
Table	19
Figures	20
Appendixes	22



Report No.: MSA-202010-409-T03

SUMMARY

The present study was to investigate the skin sensitization response of "Ear-loop" extract in guinea pigs, the testing was performed in compliance with ISO 10993-10:2010. After the treatment of the test article (polar and non-polar extraction), the extracts were applied two times for induction phase and one time for challenge phase. 24 hr and 48 hr after challenge phase, no visible changes had been found on the treated areas of either the control or the treatment group (Table 1, Appendix 1 and Figure 1), according to the criteria of "Magnusson and Kligman scale", ISO 10993-10:2010. The results indicated that the test article extract (polar or non-polar) did not cause delayed hypersensitivity on the skin of guinea pigs.



Report No.: MSA-202010-409-T03

INTRODUCTION

This study was performed in compliance with ISO 10993-10:2010 (Biological evaluation of medical devices-Part 10: Tests for irritation and skin sensitization) to evaluate the possibility of delayed hypersensitivity after topical applications with test article "Ear-loop" extract (polar or non-polar) on the skin of guinea pigs.

MATERIALS AND METHODS

- 1. Animals
 - 1.1. Species/Strain: Guinea Pig (Hartley strain)
 - 1.2. Resource: Hui Jun
 - 1.3. Body weights (Gender): 300-500 g (Male)
 - 1.4. Quarantine/acclimation (MSAT-SOP-AM-001):

Totally 30 animals were introduced in-house, animals were subjected to be quarantined and acclimated before treatment. Veterinarian ensured the animal health status before the treatment.

- 1.5. Reasons chosen for animal experimentation: Guinea pigs were proven to be suitable for skin sensitization studies.
- 2. Feeding and care (MSAT-SOP-AM-002)
 - 2.1. Housing: Guinea Pig room
 - 2.2. Environment
 - a. Temperature: 23 ± 3 °C
 - b. Humidity: $50 \pm 20\%$
 - c. Light Cycle: 12 hours light and 12 hours dark
 - 2.3. Cage and animal No.
 - a. Quarantine/acclimation: 2-3 male guinea pigs /cage
 - b. Study period: 2-3 male guinea pigs /cage



Report No.: MSA-202010-409-T03

2.4. Feed

a. Name: Lab Diet®5025

b. Brand: Lab Diet[®], U.S.A.

c. Way to supply: ad libitum

d. Source: PMI Nutrition International, U.S.A.

2.5. Drinking water

a. Sort: RO Water

b. Way to supply: ad libitum

- 3. Individual and group identification (MSAT-SOP-AM-025)
 - 3.1. Individual identification: Tested animals were identified by ear-punching.
 - 3.2. Group identification: Cages were properly labeled for identification including the Study Title/No., Administration, Observation Period, Room No., Cage No., Quantity/cage, Species, Strain, Gender, In House Date, In House Age, Animal ID No., Keeper and Deputy.
- 4. Test article and control vehicles
 - 4.1. Test article: Ear-loop
 - 4.2. Saline: CHI SHENG CHEMICAL CORPORATION, Lot. O3149;

Solutio Natrii Chloridi Isotonica 0.9% 500 ml

4.3. Cottonseed oil: SIGMA-ALDRICH, Co., Lot. MKCB9547;

CAS No.: 8001-29-4

4.4. Positive Control: 1-chloro-2, 4-dinitrobenzene (DNCB), Lot. STBG2582V

SIGMA-ALDRICH, CAS No: 101-86-0



Report No.: MSA-202010-409-T03

5. Administration of test article and control vehicles

5.1. Preparation (MAST-SOP-GE-014):

According to ISO 10993-12:2012 guideline, the test articles were extracted two times for the first and second induction phase and one time for challenge phase.

- a. Polar preparation: Measured the surface of test article, and immersed it in saline for 72±2hr at 50±2°C with constant agitation (100 rpm). The surface ratio of test article/saline was approximately 3 cm²/ml. After extraction, the extract was used immediately; the appearance of test article extract was clear and colorless without particulates present.
- b. Non-polar preparation: Measured the surface of test article, and immersed it in cottonseed oil for 72±2hr at 50±2°C with constant agitation (100 rpm). The surface ratio of test article/cottonseed oil was approximately 3 cm²/ml. After extraction, the extract was used immediately; the appearance of test article extract was clear and colorless without particulates present.

5.2. Method, route and frequency of administration:

- a. Induction phase: Two topical (upper backside) applications, intracutaneous and article spread to the dermal area. The area was from neck to scapular area.
- b. Challenge phase: One topical (lower backside) application, and article spread to the dermal area. The area was from scapular to hip area.



Report No.: MSA-202010-409-T03

5.3. Groups

Polar Group		Dose		Number of animals
"Saline"	Induction phase I	Induction phase II	Challenge phase	Male
Control-1 group	0.1 ml	0.3 ml	0.3 ml	5
Treatment-1 group	0.1 ml	0.3 ml	0.3 ml	10

Non-polar Group			Number of animals	
"Cottonseed oil"	Induction phase I	Induction phase II	Challenge phase	Male
Control-2 group	0.1 ml	0.3 ml	0.3 ml	5
Treatment-2 group	0.1 ml	0.3 ml	0.3 ml	10

6. Procedure (MSAT-SOP-ME-003)

- 6.1. Prior to the study, the furs of animal's backside were clipped from neck to scapular area with an electric animal shaver. The animals with scratches or skin diseases were removed from this study. The clipped area of left or right side was $2\times2\text{cm}^2$ each.
- 6.2. Induction phase I: On the treatment day, three kinds of solutions or emulsions were prepared from the control solution and test article extract. They were as follows:

6.2.1. Polar:

- (A)Emulsion of Freund's complete adjuvant (FCA) in saline and volume ratio of 1:1
- (B) Solution of either test article extract or Saline



Report No.: MSA-202010-409-T03

- (C)Emulsion of either test article extract or saline in E-FCA in volume ratio of 1:1
- 0.1 ml of the respective solutions/emulsions of (A), (B) and (C) of the control solution and the test article extract were injected intradermally and symmetrically into right and left scapular regions.

6.2.2. Non-Polar:

- (A)Emulsion of Freund's complete adjuvant (FCA) in cottonseed oil and volume ratio of 1:1
- (B) Solution of either test article extract or cottonseed oil
- (C)Emulsion of either test article extract or cottonseed oil in E-FCA in volume ratio of 1:1
- 0.1 ml of the respective solutions/emulsions of (A), (B) and (C) of the control solution and the test article extract were injected intradermally and symmetrically into right and left scapular regions.


Test periods	Treatment-1 group (10 animals)	Control-1 group (5 animals)
	A. Saline and FCA (1:1)	A. Saline and FCA (1:1)
Induction (I)	B. test article extract	B. Saline
	C. test article extract and E-FCA (1:1)	C. Saline and E-FCA (1:1)
Induction (Ⅱ)	test article extract	Saline
Challenge	test article extract	Saline

Test periods	Treatment-2 group (10 animals)	Control-2 group (5 animals)
	A. Cottonseed oil and FCA (1:1)	A. Cottonseed oil and FCA (1:1)
T 1 ((T)	B. test article extract	B. Cottonseed oil
Induction (I)	C. test article extract and E-FCA (1:1)	C. Cottonseed oil and E-FCA
		(1:1)
Induction (II)	test article extract	Cottonseed oil
Challenge	test article extract	Cottonseed oil




Report No.: MSA-202010-409-T03

The detail information for grouping and applied substance was tabulated above:




Polar: Treatment-1 group by 10 animals (solvent-Normal saline)




Polar: Control-1 group by 5 animals (saline)



Report No.: MSA-202010-409-T03



Non-Polar: Treatment-2 group by 10 animals (solvent-cottonseed oil)



Non-Polar: Control-2 group by 5 animals (cottonseed oil)



- Report No.: MSA-202010-409-T03
- 6.3. Induction phase II: One week later, the injection sites were applied with 10 % of sodium dodecyl sulfate (SDS) for 24 hr if there were no irritation reaction. Then, injection sites were applied by the patches that soaked with 0.3 ml of the control solution and test article extracts for another 48 hr.
- 6.4. Challenge phase: On the challenge day (two weeks after induction phase II), the furs of lower backside of the animals were clipped from scapular to hip area. An appropriate site of this hairless area was selected and applied by the patches that soaked with 0.3 ml of the control solution and test article extracts. The patches were removed 24 hr later.
- 7. Animal observations and items for examination (MSAT-SOP-ME-003)
 - 7.1. 24hr and 48hr after challenge phase, the skin response was graded according to the criteria of "Magnusson and Kligman scale, ISO 10993-10: 2010" (Appendix 2).

8. Test positive control

- 8.1. The positive control testing for sensitization assay was performed using a test sample consisting of 1-Chloro-2,4-dinitrobenzene (DNCB), which was a typically used as a positive control in this test. The DNCB concentrations used for induction and challenge in the guinea pig maximization sensitization test were 0.5% and 0.1% w/v, respectively.
- 8.2. The positive control was tested completely within three months of test article results to confirm that tested with the same source and strain of animals, as well as the same methods.
- 8.3. Test period of positive control: 10.05.2020-10.29.2020; History record from Master Lab, Report No. MSA202010.



Report No.: MSA-202010-409-T03

RESULTS

The study was to investigate the skin sensitization response of "Ear-loop" extracts in guinea pigs. Following the treatment of the test article (polar and non-polar extraction), the extracts were applied two times for induction phase and one time for challenge phase. 24 hr and 48 hr after challenge phase, no visible changes had been found on the treated areas of either the control or the treatment group (Table 1, Appendix 1 and Figure 1).

CONCLUSION

The present study was performed in compliance with ISO 10993-10: 2010 to evaluate the possibility of delayed hypersensitivity after topical applications of test article extract (polar or non-polar) on the skin of guinea pigs. After 24 hr and 48 hr challenge phase, neither the control nor the treatment group showed visible changes of skin response on the treated areas. The results indicated that the test article "Ear-loop" extract (polar or non-polar) did not cause delayed hypersensitivity in the tested guinea pigs.



Report No.: MSA-202010-409-T03

REFERENCES

- 1. Good Laboratory Practice for Nonclinical Laboratory Studies (2019) Food and Drug Administration, R.O.C.
- Good Laboratory Practice for Nonclinical Laboratory Studies. Title 21 of the U.S.
 Code of Federal Regulations, Part 58 (2019) United States Food and Drug Administration.
- 3. Biological evaluation of medical devices-Part 10: Tests for irritation and skin sensitization, ISO 10993 (2010).
- 4. Biological evaluation of medical devices-part 10: Test for irritation and skin sensitization, EN ISO 10993 (2013).
- 5. Biological evaluation of medical devices-Part 12: Sample preparation and reference materials ISO 10993 (2012).
- 6. ISO/IEC 17025:2017 General requirements for the competence of testing and calibration laboratories.
- 7. Handbook of Toxicology, Third Edition 2014, Edited by Michael J. Derelanko and Carol S. Auletta, Chapter 3. Dermal Toxicology

Report No.: MSA-202010-409-T03

Table 1. Skin Reaction in Guinea Pigs

Polar Group (Saline)	Control	Treatment
Gender	Male	Male
Number of animals	5	10
Erythema and eschar	0/5	0/10
Edema	0/5	0/10


n/n: No. of guinea pigs with abnormal clinical signs / No. of guinea pigs per group


Non-Polar Group		
(Cottonseed oil)	Control	Treatment
Gender	Male	Male
Number of animals	5	10
Erythema and eschar	0/5	0/10
Edema	0/5	0/10

n/n: No. of guinea pigs with abnormal clinical signs / No. of guinea pigs per group

Report No.: MSA-202010-409-T03


Figur1-1. Pictures for Observation of Skin Reaction (Saline)





Observation at the 24th hr after administrated Observation at the 48th hr after administrated Control-1 group

Control-1 group




Observation at the 24th hr after administrated Observation at the 48th hr after administrated


Treatment-1 group

Treatment-1 group

Figure 1-2. Pictures for Observation of Skin Reaction (Cottonseed oil)




Observation at the 24th hr after administrated



Observation at the 48th hr after administrated


Control-2 group

Control-2 goup



Observation at the 24th hr after administrated

Treatment-2 group




Observation at the 48th hr after administrated


Treatment-2 group




Report No.: MSA-202010-409-T03

Figure 1-3. Pictures for Observation of Skin Reaction Positive control group







Observation at the 24th hr after administrated Observation at the 48th hr after administrated Positive control group

History record from Master Lab, Report No. MSA202010 (Test period: 10.05.2020-10.29.2020)



Report No.: MSA-202010-409-T03

Appendix 1-1. Individual Animal Grade in Skin Response of Guinea Pigs

	Polar Groups	Gender	ID of animals	24 hr after challenge phase	48 hr after challenge phase
			1001	0	0
			1002	0	0
	Control-1 "Saline"	Male	1003	0	0
	Same		1004	0	0
			1005	0	0
			1006	0	0
			1007	0	0
			1008	0	0
	T		1009	0	0
	Treatment-1 "Ear-loop"	N.C. 1.	1010	0	0
	extract by	Male	1011	0	0
	"Saline"		1012	0	0
			1013	0	0
			1014	0	0
			1015	0	0



Report No.: MSA-202010-409-T03

Appendix 1-2. Individual Animal Grade in Skin Response of Guinea Pigs

Non-Polar Groups	Gender	ID of animals	24 hr after challenge phase	48 hr after challenge phase
0 =	, ,	1016	0	0
		1017	0	0
Control-2 "Cottonseed oil"	Male	1018	0	0
Cottonseed on		1019	0	0
		1020	0	0
		1021	0	0
		1022	0	0
		1023	0	0
Treatment-2		1024	0	0
"Ear-loop"	Male	1025	0	0
extract by	Maie	1026	0	0
"Cottonseed oil"		1027	0	0
		1028	0	0
		1029	0	0
		1030	0	0



Report No.: MSA-202010-409-T03

Appendix 1-3. Individual Animal Grade in Skin Response of Guinea Pigs

Groups	Gender	ID of animals	24 hr after challenge phase	48 hr after challenge phase
	***************************************	1006	1	1
		1007	1	_ 1
		1008	2	2
Positive control		1009	1	1
"DNCB" Induction phase	37.1	1010	2	2
(0.5% DNCB)	Male	1011	1	1
Challenge phase		1012	1	2
(0.1% DNCB)	(0.1% DNCB)	1013	1	2
		1014	1	1
		1015	2	2

History record from Master Lab, Report No. MSA202010 (Test period:10.05.2020-10.29.2020)

The positive control was tested completely within three months of test article results to confirm that was tested with the same source and strain of animals, as well as the same methods.

Appendix 2. Magnusson and Kligman Scale ISO 10993-10:2010

Patch test reaction	Grading scale
No visible change	0
Discrete or patchy erythema	1
Moderate and confluent erythema	2
Intense erythema and swelling	3



Report No.: MSA-202010-409-T03

Appendix 3. Test Article Information Sheet

Master Laboratory Co., Ltd.

Information for **☑**Test Article /□Control Article

Sponsor Company	CHIN HSIUNG FIBER CO., LTD
Sponsor Address	NO.37-12, Ta Hsin RD., Pu Yen Hsiang, Chang Hua Hsien, Taiwan
Contract study item	☑ Base on the contract ☐ Others:
Name of test article	Ear-loop
Major components	Spandex Yarn+ Nylon Yarn
Sample status	□Sterilized (□Gamma □EO □Steam) ☑ Not Sterilized
Storage condition	✓ Room temperature (10°C~30°C)
Expiry day	
Specific requirement	
Batch/ Lot number	 ☑ Base on the specific number on the package: ☐ Base on the date on the package ☐ Base on the arrived date ☐ Others:
Extract by	Weight (0.2g/ml) Total weight of each test article:
	Surface (Sample thickness: 🗹>1.0mm 🔲 0.5-1.0mm 🖂<0.5mm) Total area surface of each test article:
Absorption	✓ Non absorption
	☐ Water absorption rate: / Oil absorption rate:
Sponsor Signature	Hank Shih



Report No.: MSA-202010-409-T03

Appendix 4. Test Article

